



# More Than Mobility Inc

5282 Catalina Drive, Nanaimo, BC, V9V 1H1 Email: [morethanmobility@shaw.ca](mailto:morethanmobility@shaw.ca)

Tel Toll Free 1-866-339-0081 or (250) 668-0081



## 4LA 202 Exoskeletal LA Ankle

The Exoskeletal LA Ankle is designed for inclusion with the final lamination.

The ribbed surface is encased with the lamination.

Diameter – 42 mm

Height - 118 mm

Weight – 393 grams