

LA Ankle 201



The LA Ankle is intended for water leisure activity and is primarily a swimming aid.

- Easy locking system
- Once attached the ankle closely resembles a natural ankle joint.
- It can be locked into 70° of plantar flexion for swimming

It is not intended as a functional ankle for day to day walking. The user should use their normal prosthesis if intending to do walking other than transferring to and from the pool, lake or ocean.

Tested to weight range 100kg/ 220 Lbs

- Fitted with a titanium pyramid providing $\pm 4^\circ$ angulation and 180° rotation.
- Stainless Steel Lock Pin
- Main body is Aluminum, hard-coated black to prevent corrosion.



Weight:

Unit weighs **only** 410 grams / 14 Oz

Warranty:

The unit has a two warranty.



More Than Mobility Inc

More Than Mobility Inc. 5282 Catalina Drive, Nanaimo, British Columbia, Canada. V9V 1H1

1-866-339-0081 or (250) 668-0081 morethanmobility@shaw.ca www.morethanmobility.ca