

The Pneumatic Post-Amputation Mobility Aid

Enhancing early patient rehabilitation

The Pneumatic Post-Amputation Mobility Aid (PPAM Aid) is a device designed for partial weight bearing rehabilitation. Early use of the PPAM Aid assists in the reduction of residual limb oedema and promotes mobility, offering clear psychological benefits for the patient.

The PPAM Aid kit consists of several components:

- Frame
- Non-tacky stump and cushion bags
- Pump
- Quick release valves for bags
- Shoulder support straps
- Crucible sling support



The device is suitable for below knee, knee disarticulation and long above-the-knee amputees.

Benefits

- Prepares residual limb for wearing a prosthesis
- Allows for the re-education of postural reactions, balance and gait
- Aids reduction of residual limb oedema
- Range of frame sizes able to cater for bariatric patients
- Fast-acting pump
- Can be used over soft dressing, bandaging or a plaster cast

Ordering information

PPAM Aid kit part numbers

Part number	Description
PPAM-REHABKIT	PPAM Aid rehabilitation kit (includes 6 x frames, 6 x bags, pump and accessories)
PPAM-BARKIT	PPAM Aid bariatric kit (includes 3 x frames, 3 x bags, pump and accessories)

PPAM Aid componentry part numbers

Part number	Description	Latex free
Frames		
6164	650mm frame	
6163	750mm frame	
6154	850mm frame	
6164/L	650mm frame (Large)	
6163/L	750mm frame (Large)	
6154/L	850mm frame (Large)	
Bags		
6174	Below Knee stump bag	X
6174-N	Non-tacky Below Knee stump bag	X
6175	Above Knee stump bag	X
6175-N	Non-tacky Above Knee stump bag	X
6148	Cushion bag	X
6148-N	Non-tacky cushion bag	X
6174/L	Below Knee stump bag (Large)	X
6174/L-N	Non-tacky Below Knee stump bag (Large)	X
6175/L	Above Knee stump bag (Large)	X
6175/L-N	Non-tacky Above Knee stump bag (Large)	X
6148/L	Cushion bag (Large)	X
6148/L-N	Non-tacky cushion bag (Large)	X
Pump		
6253	PPAM Aid pump set	
4626	Pump rubber bulb replacement (for use with old pump, part number 6153)	
Accessories		
4599	650mm to 750mm crucible sling support	
4630	850mm crucible sling support	
6178-9	Shoulder support straps	
4603/L	Replacement rubber rocker	
6181	Quick release male valve for existing PPAM bags	
6181-F	Quick release female valve for existing pump solutions	
6181-KIT	Quick release valve kit (includes 1 x 6181-F and 6 x 6181)	